

## Northern blotting

**Reagents:** Agarose GTG: Cat # 50070, FMC BioProducts (Rockland, Maine)  
10xMOPS:  
37% formaldehyde  
DEPC water  
Ethidium bromide  
Deionized formamide

### Procedure:

#### 1. Gel

horizontal gel 14 slots, wide comb, 3.5ug per lane

Agarose gel	agarose GTG	1.2g
	H <sub>2</sub> O	87ml in Erlenmeyer
	Microwave until dissolved (~2min)	
	10xMOPS	10ml
	cool to touch	
	37% formaldehyde	5.1ml

Mix above solution gently to avoid forming air bubbles. Pour all into gel box, wait until it cools to form a 0.5cm gel

#### 2. Samples

a. Making buffer:	Sample buffer	a. formamide (deionized)
500ul		
	10x MOPS	150ul
	DEPC water	150ul
	40% formaldehyde	200ul
	b. "northern dye"	200ul

Mix A and B, store at -20°C

Diluted ethidium bromide	ethidium bromide	5ul
	DEPC water	1ml

#### b. Loading samples

Label the sample tubes. Add 12ul sample buffer and 3ul diluted ethidium bromide (5 ul of 10 mg/ml in 1000 ul water) to all the tubes. Add sample and DEPC water according to below table. Heat the sample tubes at 65°C for 5 minutes. Load all the sample (~25ul) into sample well (see map show bellow). Run the gel at 50Volts for 4-5 hours.



### 3. Blotting

- a. Rinse the gel twice with 10xSSC for 15 min. to remove formaldehyde, which increases transfer efficiency of the gel).
- b. Put blotting pad in a plastic container. Wet top 2-3 blotting pads with 10xSSC.
- c. Cut Zeta-Probe GT membrane to gel measurements (11x14 inches), wet with 10xSSC.
- d. Place the gel "upside-down" (RNA sink on the bottom of the gel) on wet pads, lay membrane on it, remove air bubbles between membrane and gel, put 2 wet pads above and cover with dry pads.
- e. Transfer for overnight.
- f. Observe the transfer result under UV. The bands are easily observed on the reverse side of RNA on the membrane. For a good transfer, no bands will be left in the gel.
- g. Link the RNA to membrane (stratalink) in UV-oven (autolink+start).
- h. Put the membrane into a pouch, store it at 4°C.

### 4. Prehybridization

- a. Hybridization solution Hybrisol (ONCOR, Gaithersburg, MD)
- b. Preheat hybrisol at 55°C (dissolve all ingredients) until it is clear.
- c. Put membrane into hybridization tube and add 15ml hybrisol.
- d. Prehybridization at 41°C in oven with rotation for at least 3 hours.

5. Radiolabel probe
- |   |            |
|---|------------|
| DNA   | 25ng (1ul) |
| Water   | 22ul       |
| Random Oligonucleotide primers  | 10ul       |
| Boil for 5 minutes at 95-100°C, spin-down the liquid                                |            |
| 5x Primer buffer  | 10ul       |
| □CTP <sup>32</sup>  | 5ul        |
| T7 DNA polymerase   | 1ul        |
| Incubate for 5 minutes at 37°C  |            |
| Blue sample buffer (optional)<br>to stop the reaction<br>(any labeling kit will do) | 10ul       |

### 6. Purify the probe

STE Select-D, G-50 spin column

- a. Invert column several time to fully resuspend gel, Remove top closure first, and then bottom one

- b. Let the buffer drain, put it in a collection tube, and spin at 1000rpm for 1.5 min.
- c. Change the collection tube, and discard the old one.
- d. Put all the sample (~ 60ul) onto the center of gel bed.
- e. Spin at 1000rpm for 1.5 min.
- f. Check probe with Geigercounter: put 1ul of 4 ml of scintillation solution and count cpm. Use about  $1 \times 10^7$ /blot

## 7. Hybridization

- a. Put radioactive probe into a 0.5ml screw cap tube, boil the probe ( $1-2 \times 10^6$  cpm/ml Hybrosol) for 3 minutes in a heating block, put it into the hybridization tube with membrane, and put that into the hybridization oven.
- b. Incubate/rotate at  $42^\circ\text{C}$  for overnight.
- c. Wash the membrane 2x for 15 minutes in 2x SSC/0.1% SDS at room temperature. This can be done in the hybridization oven.
- d. Wash the membrane 2x for 20 minutes in 0.1x SSC/0.1% SDS at  $68^\circ\text{C}$ . This is also done in the hybridization oven

## 8. Autoradiography//other imaging analysis

- a. Let membrane dry briefly, wrap it, and put it in an X-ray cassette with and intensifying screen.
- b. Incubate at  $-70^\circ\text{C}$  (time depends on how “hot” the signal is by gieger). Develop the X-ray film in the developer on the 1<sup>st</sup> floor.
- c. Alternatively, place the membrane (covered with plastic wrap) in molecular dynamics imaging system for 1 h to 1 day.