

## 1. RT-PCR

### 2. RT with random primers

- a. Mix and briefly centrifuge each component before use
- b. Add 2ul (1-5ul) sample RNA or the same volume of DEPC H<sub>2</sub>O to a sterile tube.
- c. Prepare a master mixture in sterile 0.5ml tubes as follows
- d.

Total sample

1ul 10 mM dNTPs

1ul Random primer

6ul DEPC treated H<sub>2</sub>O

- e. Add 8 ul the master mixture to each sample tube, incubate at 65<sup>0</sup>C for 5 minutes, and then place on ice for at least 1 min.
- f. Prepare the following reaction mixture (adding each component in the indicated order)

2ul 10x RT buffer

4ul 25 mM MgCl<sub>2</sub>

2ul 0.1 M DTT

1ul Rnasin (or similar RNase inhibitor)

Add 9ul of reaction mixture to each RNA/primer mixture, mix gently, and briefly centrifuge.

- g. Incubate at 25<sup>0</sup>C for 2 min.
- h. Add 1ul of SUPERScript II RT to each tube, mix and incubate at 25<sup>0</sup>C for 10 min.
- i. Transfer the tubes to 42<sup>0</sup>C and incubate for 50 min.
- j. Terminate the reaction at 70<sup>0</sup>C for 15 min and chill on ice.
- k. Briefly centrifuge the sample down to the bottom and add 1ul of Rnase H to each tube and incubate for 20 min at 37<sup>0</sup>C.
- l. Add 30ul DEPC H<sub>2</sub>O, and store at - 20<sup>0</sup>C before PCR. (optional based on expression level of product)

## 3. PCR

- a. Prepare the following reaction mixture

5ul cDNA

5ul 10xPCR buffer 15mM MgCl<sub>2</sub> (Perkin Elmer)

10ul 1 mM dNTPs

1 ul 0.1ug/ul sense primer  
1 ul 0.1ug/ul anti-sense  
27.5 ul distilled water

- b. Add 44.5ul PCR mix to each tube.
- c. Add 5ul cDNA samples and a drop of mineral oil.
- d. Put the sample tubes into PCR machine.
- e. Start file # 62 for a hot start at 95<sup>0</sup>C for 10 minutes.
- f. Stop file #62, add 0.5ul of AmpliTaq Gold polymerase (Perkin Elmer).
- g. Enter file #75 and start.