

Transformation of Plasmids:

If DNA is on a piece of filter paper cut it out and put into a test tube with 100 ul water/TE (10mM tris, 0.1mM EDTA)

Locate or buy some competent E-coli. We use Nova blue (In-vitrogen), XL-blue (stragagene form the stockstore), or topo (from TA cloning kit, also from In-vitrogen). They are stored at –80 in 20 ul aliquots

Add 1 ul of DNA to one tube of E-coli. Controls are: no DNA or control plasmid DNA (like PUC, Bluescript or similar)

Incubate cells and DNA 30 minutes on ice.

Heat shock ~45-90 seconds at 41 C

Add 100 ul of “SOC” (from stock, a yellow nutrient media) and shake at 37 C (to recover)

Plate 10 ul and 90 ul separately in LB (luria broth) agar plates containing 50ug/ml ampicillin and incubate the plate upside down 37C overnight

Select colonies for amplification in LB broth containing amp (stock is located –20C, dilute 1/1000, ie for 500 ml media use 500 ul).

Extract DNA from overnight cultures using a mini (5 ml) midi (50 ml) or maxi (500ml) plasmid preparation kit. Before spinning cells (as per the protocol in the kit) save 700 ul and add 300 ul 50% glycerol in a nunc tube. Save this stock at –80 for future use.

Optional: squirt Xgal (50 ul of 10% in n-n dimethlyformamide to dissolve) and 50 ul IPTG (0.1 M in water) for “blue white selection” onto the plates and gently rub in with a flamed bent glass rod.

Determine if the DNA you transformed is correct by:

- 1) Studying the restriction map of the clone and select diagnostic restriction enzymes to test. Run a 1% agarose gel and stain with ethidium bromide.**
- 2) PCR the DNA using universal primers that flank the insertion site (multipurpose cloning site) using restriction enzyme digest.**
- 3) Send plasmid to Nikki for DNA sequencing using primers adjacent to the cloning site and compare the returned sequence by “blast” analysis (<http://www.ncbi.nlm.nih.gov/gorf/bl2.html>).**